


McDowell + Service
Dental Laboratory Ltd.



Crown & Bridge Comparison Chart

	Standard PFM	Premium PFM	Full Metal Crown	Fortress Zirconia	Full Zirconia	Chameleon Zirconia	e-MAX
Material Used	Non-precious coping with porcelain overlay	Non-precious, precious alloy with layered porcelain	Non-precious or precious alloy	Solid Zirconia restoration	Zirconia coping with multi-layered porcelain overlay	Press lithium disilicate stained or multi-layered porcelain overlay	
Method	Cast or laser sintered and layered porcelain	Cast or laser sintered and layered porcelain	Cast for precious or laser-sintered	Milled using CAD/CAM technology	Designed & milled using CAD/CAM technology and built up by hand	Pressed or designed & milled using CAD/CAM technology & built up by hand	
Aesthetics	☆☆☆	☆☆☆☆	☆	☆☆☆	☆☆☆☆☆☆	☆☆☆☆☆☆	
Strength	☆☆☆☆	☆☆☆☆	☆☆☆☆☆☆	☆☆☆☆☆☆	☆☆☆☆	☆☆☆☆	
Add-On Options (Fit to denture, attachments, etc. Charges may apply)	✓	✓	✓	✓	✓	✓	
Benefits	Traditional method	Traditional method, choice of alloys & multi-layer build	High durability, strong, precious/non precious alloy	Good aesthetics, high durability, extreme strength, metal-free & bio compatible	Premium aesthetics, excellent fit, metal-free & bio compatible	Premium aesthetics, excellent fit, metal-free & bio compatible	
Recommended usage	Anterior and posterior	Anterior and posterior	Posterior teeth	Posterior teeth & Bruxers	Anterior and posterior	Anterior and posterior	
Recommended bonding	Traditional cementation or bonding technique	Traditional cementation or bonding technique	Traditional cementation or bonding technique	Traditional cementation or bonding technique	Traditional cementation or bonding technique	Bonding technique ONLY	
Warranty	✓	✓	✓	✓	✓	✓	

McDowell + Service
Dental Laboratory Ltd.

4c Dargan Court Dargan Crescent Belfast BT3 9JP

Tel. 02890 370 708

Email: enquiries@mcdowell-service.com www.mcdowell-service.com

