



# VITA MFT

Multi Functional Teeth

Mould chart acrylic teeth




Date of issue 01.13


VITA shade, VITA made.

**VITA**

Available in VITA SYSTEM 3D-MASTER  
and VITA classical A–D shades



Ovoid



Triangular

**040**



**T41**



**044**



**T43**




**049**



**T46**



**T49**



**VITA MFT – The new dimension in aesthetics**

The natural-looking anatomy of VITA MFT anteriors, with approximal contacts in the area of the incisal third being slightly placed towards the palate, enables a trouble-free and appearingly organic alignment of the teeth.

The anatomically correct labial curve leads to a physiologically perfect set-up of the teeth where discretely broadened tooth necks support an age-appropriate papilla design.

## Lower anteriors



Rectangular

**R42**



**R45**



**R47**



Square

**S47**



**S50\***



**L33**



**L34**



**L37**



**L37LN**



**L40**



The mould **L37LN** (Long Neck) was specifically developed for simulating particular periodontal characteristics which, for example, occur in partial prosthetics.

**PU29 – PL29**



**PU31 – PL31**



**PU33 – PL33**



**easy-centric® – the new, effortless way for setting up dentures**

The occlusal surface design of VITA MFT is based on the innovative principle called easy-centric which supports a professional and economical set-up procedure. Moreover, it leads to increased productivity in a quality well adapted to the patient's needs.

Nr.		Nr.	mm	
	mm		Nr.	mm
<b>O40</b>	40.0	<b>L33</b> 32.6	<b>PU29</b> 28.8	<b>PL29</b> 29.6
<b>O44</b>	43.9	<b>L34</b> 34.0	<b>PU29</b> 28.8	<b>PL29</b> 29.6
<b>O49</b>	49.2	<b>L40</b> 39.5	<b>PU33</b> 33.0	<b>PL33</b> 33.6
<b>T41</b>	40.5	<b>L33</b> 32.6	<b>PU29</b> 28.8	<b>PL29</b> 29.6
<b>T43</b>	43.3	<b>L34</b> 34.0	<b>PU31</b> 30.6	<b>PL31</b> 31.1
<b>T46</b>	46.1	<b>L37</b> 36.8	<b>PU31</b> 30.6	<b>PL31</b> 31.1
<b>T49</b>	48.8	<b>L40</b> 39.5	<b>PU33</b> 33.0	<b>PL33</b> 33.6
<b>R42</b>	42.1	<b>L34</b> 34.0	<b>PU29</b> 28.8	<b>PL29</b> 29.6
<b>R45</b>	44.5	<b>L37</b> 36.8	<b>PU31</b> 30.6	<b>PL31</b> 31.1
<b>R47</b>	46.6	<b>L37</b> 36.8	<b>PU33</b> 33.0	<b>PL33</b> 33.6
<b>S47</b>	46.8	<b>L40</b> 39.5	<b>PU33</b> 33.0	<b>PL33</b> 33.6
<b>S50</b>	49.7	<b>L40</b> 39.5	<b>PU33</b> 33.0	<b>PL33</b> 33.6



### VITA MFT – Anteriors and posteriors

The whole range of indications in the pleasantly small VITA MFT product line. A reduced and harmonious range of shades/moulds ensures for increased economic efficiency.

The multifunctional VITA MFT tooth is suited for:

- Full dentures
- Implant prosthetics
- Telescope and cone constructions
- Attachment Dentures
- Model cast prostheses
- Interim solutions

### The shade range of VITA MFT covers

- 11 VITA SYSTEM 3D-MASTER shades


VITA SYSTEM 3D-MASTER										
1M1	2L1.5	2M1	2M2	3L1.5	3L2.5	3M2	3R2.5	4L1.5	4M2	5M1
~B1	~B2	~A2	~C2 ~D4	~B3	~A3.5	~A4	~C4			

- 13 VITA classical A–D shades

VITA classical A–D												
A1	~A2	A3	~A3.5	~A4	~B1	~B2	~B3	~C2	C3	~C4	D3	~D4
	2M2		3R2.5	4M2	1M1	2L1.5	3L2.5	3L1.5		5M1		3L1.5



With the unique VITA SYSTEM 3D-MASTER, all natural tooth shades can be systematically determined and perfectly reproduced.



Please note: Our products must be used in accordance with the instructions for use. We accept no liability for any damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of application. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers that are not compatible or not authorized for use with our product. Furthermore, our liability for the accuracy of this information is independent of the legal basis and, in as far as legally permissible, shall always be limited to the value as invoiced of the goods supplied, excluding value-added tax. In particular, as far as legally permissible, we do not assume any liability for loss of earnings, indirect damages, ensuing damages or for third-party claims against the purchaser. Claims for damages based on fault liability (culpa in contrahendo, breach of contract, unlawful acts, etc.) can only be made in the case of intent or gross negligence. The VITA Modulbox is not necessarily a component of the product.

Date of issue of these directions for use: 01.13

VITA Zahnfabrik is certified according to the Medical Device Directive; the following products bear the mark  0124 : **VITAMFT**

US 5498157 A · AU 659964 B2 · EP 0591958 B1

# VITA

VITA Zahnfabrik H. Rauter GmbH & Co. KG  
Postfach 1338 · D-79704 Bad Säckingen · Germany  
Tel. +49(0)7761/562-0 · Fax +49(0)7761/562-299  
Hotline: Tel. +49(0)7761/562-222 · Fax +49(0)7761/562-446  
[www.vita-zahnfabrik.com](http://www.vita-zahnfabrik.com) · [info@vita-zahnfabrik.com](mailto:info@vita-zahnfabrik.com)  
 [facebook.com/vita.zahnfabrik](https://facebook.com/vita.zahnfabrik)