

## McDowell + Service Guide to All Ceramic Restorations

Restoration Type	Restoration Material				
	Resin Bonded	PFZ	e.Max	e.Max Zirpress	Fortress Zi
Crown	Yes	Yes	Yes	Yes	Yes
Bridge	No	Yes	Up to 3 units	Up to 6 units	Yes
Inlay or Onlay	Yes	No	Yes	No	Possible
Maryland	No	Possible	No	Possible	Possible
Cantilever Bridge	No	Yes	No	Yes	Yes
Inlay or Bridge	Possible	Possible	Possible	Yes	Possible
Veneer	Yes	No	Yes	No	No
Implant Abutment	No	Yes	Possible	No	No
Implant Bridge Screw Retained	No	Yes	Possible	Up to 6 units	Yes
Implant Bridge Cement Retained	No	Yes	Up to 3 units	Up to 6 units	Yes
Implant Crown Screw Retained	No	Yes	Possible	Yes	Yes
Implant Crown Cement Retained	No	Yes	Yes	Yes	Yes