

# natura



The Quality tooth at a mid-range price

**MOULD CHART**

Upper Anterior Teeth

**K** TAPERED



**S** SQUARE



**K21** | | | | 40.2 | 9.4 | 7.5

**S11** | | | | 41.5 | 8.6 | 8.1

**K22** | | | | 43.0 | 9.6 | 7.9

**S12** | | | | 45.1 | 8.5 | 8.6

**K24B** | | | | 43.8 | 10.6 | 8.5

**S13** | | | | 43.1 | 9.4 | 8.4

**K25** | | | | 47.4 | 10.5 | 8.9

**S14** | | | | 46.6 | 10.7 | 8.9

**K27** | | | | 49.2 | 10.9 | 9.1

**S15** | | | | 50.3 | 10.7 | 9.4

**K41** | | | | 42.7 | 9.8 | 8.2

**S17** | | | | 52.0 | 12.8 | 9.8

**K42** | | | | 44.3 | 9.6 | 8.0








**S66** | | | | 46.4 | 10.3 | 8.8

# Lower Anterior Teeth

## O OVAL








## LOWER



## Posterior Teeth

## Combination Tables

### POSTERIOR



Anterior Teeth		Posterior Teeth
Upper	Lower	
S11	L3	P3
S12	L5	P3
S13	L5	P3/P5
S14	L6	P5
S15	L8	P4
S17	L9	P6
S66	L5	P5
K21	L3	P3
K22	L3	P3/P5
K24B	L4	P2
K25	L8	P4
K27	L8	P4
K41	L3	P3/P5
K42	L5	P3/P5
O31	L3	P3
O32	L3	P3/P5
O34	L5	P5
O36	L7	P5
O54	L5	P5

Instructions for obtaining optimum results are available from us on request.



# natura

the natural choice

After more than five years research and development, Schottlander have created **Natura** - a range which is vital and natural in appearance yet priced as a mid range tooth.

**Natura** anterior teeth use skillful shading, internal mamelon effects and a subtle response to changing light conditions to achieve their remarkably lifelike appearance.

**Natura** posterior teeth - with natural translucent cusps - are manufactured to the same high aesthetic standard and are both more wear resistant than standard PMMA teeth and exceptionally strong and durable in the mouth.



#### Their many advantages include:

- Internal mamelon effects in both upper and lower centrals for a more natural appearance.
- Naturally occurring blue/white opalescence and fluorescence especially visible in upper laterals.
- Anterior aesthetics carried through to posterior teeth.
- CE marked and manufactured under ISO 9001 quality standard and in accordance with EN ISO 22112 Standard for Artificial Teeth.
- Excellent bond to acrylic denture base material.
- High solvent and heat resistance.

**Natura** the quality tooth at a mid-range price.



In 2004 Schottlander were honoured to receive the Queen's Award for the Enigma and Natura teeth and denture systems.

## Schottlander's Innovation Recognised

The official citation from the Queen's Awards Office read:

"The award to Davis Schottlander & Davis Ltd recognises this company's success in continuously improving a range of teeth and dentures.

The Company has created Enigma and Natura, two ranges of teeth and dentures that not only look more natural but also function better. This has been achieved by using advanced acrylic materials, a wider range of colours and better mechanical design, to give an improved "bite". Options available now take account of the different characteristics of ethnic groups. Continuous innovation in design and professional marketing have enabled this award winner to achieve a substantial share of the market and considerable commercial success whilst providing customers with better looking and performing dentures".

Tel +44 (0)1462 480848 Fax +44 (0)1462 482802  
Freephone 0800 97 000 79 Freefax 0800 97 000 78  
e-mail: [sales@schottlander.co.uk](mailto:sales@schottlander.co.uk)  
[export@schottlander.co.uk](mailto:export@schottlander.co.uk)

***schottlander***  
*committed to quality*

Davis Schottlander & Davis Ltd  
Fifth Avenue,  
Letchworth Garden City  
Herts, SG6 2WD England