

# BIOCOMPATIBILITY CERTIFICATE

**Testmaterial:** REMATITAN® cast metal  
Titan grade 1 according to DIN 17850/ISO 5832-2/ASTM F67

**Supplier:** DENTAURUM J.P. Winkelstroeter KG  
Turnstraße 31 D-75228 Ispringen

**Study performed:**

The following study was performed in order to determine the biocompatibility of the dental material according to ISO 10993-5, 1992: „Test on cytotoxicity; *in vitro* methods“. The material was produced at DENTAURUM.

### CYTOTOXICITY

The cytotoxic effect of the test material was examined in a cell culture test *in vitro* using L 929-cells (fibroblast cell line). Examination of materials for medical devices with the extraction method and protein determination.

**Result:** REMATITAN® did not show any cytotoxic effect.

**anawa BIOSERVICE Scientific Laboratories GmbH Munich**


Behringstraße 6

D-82152 Planegg



Dr. Albrecht Poth

Date: April 26, 1996



Dr. Naveed Honarvar

Date: April 26, 1996